

# FUSSEN by orangedental

## Intraoral scanner model S6500

### Model S6500

NEW



- Large field of view 16 x 13 mm
- Complete jaw scan in less than 2 minutes
- High accuracy only 10 µm tolerance
- Super light and compact only 220 g
- 3 standard and 1 mini scan tip
- New ScanFlex software with many AI functions



- New, modern and stable cart with mini PC, touch screen and emergency power supply



# FUSSEN by orangedental

## Intraoral scanner

### Model S6500

	S6500
Type	Cable
Scanprinciple	Triangulation
Light source	Structured light
Complete jaw scan (min)	< 2
Scanning speed (FPS)	20-25
Depth of field (mm)	0-15
Field of view (FOV)	16x13/12x9
Anti-fogging (min)	< 1
Weight (g)	220
Size (mm)	240x45x36
G-Motion	Yes
Precision single crown	10µm (±1,1µm)
Size of the scanning tips (mm)	Standard: 95,3 x22,5x16,6 Mini: 94x 17,9x 13,2
Included scan tips	3 standard scanning tips 1 Mini scanning tip
RRP plus VAT.	12.900 €

